

Pharmaceutical Solution, Inc.

High Shear Granulator

Application:

Completely self contained bottom drive high shear granulator.

Features:

- 1. Lab and Pilot scale models have interchangeable product bowls to accommodate different batch size.
- Three-level security authorization, parameters controlled by PLC, auto-save for all data which can be exported by USB devices.
- 3. Liquid Adding: Airless atomization, wetting agents are aerosolized airlessly, the droplet becomes more uniform, smaller and even. This can effectively solve highly viscous materials
- 4. End Point Determination: Equipment can be customized, customizable end point output of torque, power, electric current, etc.
- 5. Food grade skeleton oil sealed rotating shaft design; it eliminates dust bonding.
- 6. The gap between the impeller and the bottom of the bowl is minimized (0.3-1mm), which can solve the problem of material deposited in the bottom of the bowl.
- 7. Mixing bowl, impeller and chopper parts can be changed in 3 minutes without any tools.

Address: 46500 Fremont Blvd, STE 702, Fremont, CA94538

Phone Number: (510) 573-0562







High Shear Granulator

Lab to Pilot Scale Models and Technical Data

| | Mini G | G6 | G10 | G20 | G30 | | | | |
|--------------------------------|-------------------------------|----------------------------|------------------------|------------------------|----------------------------|--|--|--|--|
| | (Interchangeable bowl) | (Interchangeabl e bowl) | (Interchangeable bowl) | (Interchangeable bowl) | (Interchangeabl e bowl) | | | | |
| Bowl Volume (L) | 0.1、0.2、0.5 | 0.5、1、3、6 | 1、3、6、10 | 1、3、6、10、 20 | 6、10、20、30 | | | | |
| Working Capacity (kg/batch) | 0.02-0.10 | 0.05-2.0 | 0.1-3.0 | 0.1-6.0 | 2-10 | | | | |
| Impeller Speed (rpm) | 50-1400 | 50-480 | 50-480 | 50-480 | 50-480 | | | | |
| Chopper Speed (rpm) | 300-3000 | 300-3000 | 300-3000 | 300-3000 | 300-3000 | | | | |
| Power | 220V/ 1 phase 3 wire /50Hz | 380V/ 3 phase 5 wire/50Hz | | | | | | | |
| Total Power (kw) | 0.36 | 2 | 3 | 4.5 | 4.5 | | | | |
| Air Pressure (MPa) | n/a | · 0.3 | → 0.3 | · 0.3 | · 0.3 | | | | |
| Dimensions (mm) | 270×300×300 | 1200×550×1100 | 1350×600×1100 | 1200×600×1150 | 1200×800×1200 | | | | |

Pilot to Production Scale Models and Technical Data

| | G50 | G100 | G200 | G300 | G600 | | | |
|-----------------------------|---------------------------|------------------------|--------------------|--------------------|--------------------|--|--|--|
| | (Interchangeable bowl) | (Interchangeable bowl) | (fixed bowl) | (fixed bowl) | (fixed bowl) | | | |
| Bowl Volume (L) | 3、6、10、20、 50 | 20、50、100 | 200 | 300 | 600 | | | |
| Working Capacity (kg/batch) | 0.1-20 | 10-40 | 20-80 | 30-120 | 60-240 | | | |
| Impeller Speed (rpm) | 50-480 | 50-480 | 50-480 | 30-300 | 30-300 | | | |
| Chopper Speed (rpm) | 300-3000 | 300-3000 | 300-3000 | 300-3000 | 300-3000 | | | |
| Power | 380V/3 phase 5 wire/ 50Hz | | | | | | | |
| Total Power (kw) | 7.5 | 16 | 21 | 30 | 42 | | | |
| Air Pressure (MPa) | > 0.4 | · 0.5 | · 0.5 | · 0.5 | · 0.3 | | | |
| Dimensions (mm) | 1800×800×1300 | 2100×1200×1700 | 2500×1500×1 800 | 2800×1800×2 000 | 2800×1800× 2000 | | | |

Address: 46500 Fremont Blvd, STE 702, Fremont, CA94538

Phone Number: (510) 573-0562

EPS High Shear Granulator Provides:

- Batch Size Flexibility
- Easy Scalability
- Optimum Process Control
- Excellent Handling and Friendly User Technology
- Easy Changeover and Cleaning